

MOBIDERM® Intimate Bra



☐ Order ☐ Quotation ☐ Renewal


Date: ____ / ____ / ____

Patient Name: _____

Medical Professional Name: _____

Fitter Name: _____ Case No. for renewal: _____


☐ 1st Treatment Phone: _____ Email: _____



1. Band Size

Measure the circumference just under the breast for band size.

_____ (cm)



2. Cup Size


Measure the circumference of the fullest part of the breast for cup size.

_____ (cm)

Sizing Chart		XS	S	M	L	XL
① Band Size	CM	68 - 77 cm	78 - 87 cm	88 - 97 cm	98 - 107 cm	108 - 117 cm
	IN	26.8 - 30.3 in	30.7 - 34.3 in	34.6 - 38.2 in	38.6 - 42.1 in	42.5 - 46 in
② Cup Size	AA/A	80 - 89 cm 31.5 - 35 in	90 - 99 cm 35.4 - 39 in	100 - 109 cm 39.4 - 43 in	112 - 119 cm 44.1 - 46.9 in	120 - 129 cm 47.2 - 50.8 in
	ITEM #	☐ T571301201AAA24	☐ T571301202AAA24	☐ T571301203AAA24	☐ T571301204AAA24	☐ T571301205AAA24
	B/C	84 - 93 cm 33.1 - 36.6 in	94 - 103 cm 37 - 40.6 in	104 - 113 cm 41 - 44.5 in	116 - 123 cm 31.5 - 45.7 in	124 - 133 cm 48.8 - 52.4 in
	ITEM #	☐ T5713012010BC24	☐ T5713012020BC24	☐ T5713012030BC24	☐ T5713012040BC24	☐ T5713012050BC24
	D/E	88 - 97 cm 34.6 - 38.2 in	98 - 107 cm 38.6 - 42.1 in	108 - 117 cm 42.5 - 46.1 in	120 - 127 cm 47.2 - 50 in	128 - 137 cm 50.4 - 54 in
	ITEM #	☐ T5713012010DE24	☐ T5713012020DE24	☐ T5713012030DE24	☐ T5713012040DE24	☐ T5713012050DE24
HCPCS: A6528		Description: Gradient compression garment, bra, for nighttime use, each			Qty: _____ (includes 2 Mobiderm Breast Pads)	

*Rounding may apply to inch measurements

Additional **Pre-made Mobiderm® Breast Pads** available for order. The pads accommodate all cup sizes and come in Small Block only.



	Item Number	Size	X	Y	Z
☐	T57135120199970	XS	260 mm	155 mm	100 mm
☐	T57135120299970	S	285 mm	160 mm	100 mm
☐	T57135120399970	M	310 mm	165 mm	105 mm
☐	T57135120499970	L	335 mm	170 mm	110 mm
☐	T57135120599970	XL	360 mm	180 mm	120 mm
Suggested HCPCS *: A6609, A6593, or A6519				Qty: _____	

Additional **Mobiderm® Mobilizing Pads** available in Big and Small block. Sold separately to cut and use for different desired shapes or sizes.

Comments: _____

*The Suggested HCPCS varies based on usage with bandaging (A6609), daytime garments (A6593), or nighttime garments (A6519).

	Pad Size	Block Size	Product Code
☐	25cm x 1m	Small Block 5mm x 5mm	37100103310001
☐	20cm x 1m	Big Block 15mm x 15mm	37100103310002
☐	25cm x 25cm	Small Block 5mm x 5mm	37100102502501
Suggested HCPCS *: A6609, A6593, or A6519			Qty: _____